

Permobil power standing

Evaluation checklist

A qualified medical professional **must** determine and document in the medical record that the individual is safe and medically cleared to stand.

Equipment

Permobil F5 Corpus VS or M Corpus VS	Metric Allen keys
Permobil QuickConfig app (Power Platform equipped wheelchair)	Single switches
QuickConfig quick start guide (Power Platform equipped wheelchair)	Lap belt and chest bar/strap
Permobil programming key and ICS switchbox (non Power Platform equipped wheelchair)	Anterior knee support
Permobil power standing programming guide (non Power Platform equipped wheelchair)	Cushion & back support that are appropriate for client
Rigid measuring tool	Arm support appropriate for client
	Drive control appropriate for client

Instructions

Step 1: Take chair through stand sequence before transferring the client and ensure it is set through checkpoint programming within clinician-evaluated ROM tolerances

Step 2: Determine how the client or caregiver will control the standing sequence (joystick, alt drive, single switch, ICS)

Step 3: Prior to placing client in the chair, adjust seat depth to client's measurements based on clinician's evaluation. Confirm programming for seat depth matches the actual, adjusted seat depth

Step 4: Lower leg length: adjust footplate height to client's measurements based on clinician evaluation

Step 5: Adjust footplate angle to accommodate seated ankle range based on clinical evaluation

Step 6: Permobil chest support is required for standing. Height, depth, and angle is determined by clinical evaluation and stability needs

Step 7: Permobil anterior knee support is required for standing. Single post anterior knee support depth should have 1"/2.5 cm of space when in a seated position. The swing-away anterior knee support should make loose contact with the body when in the seated position

Step 8: Additional positioning components as determined by evaluation: ex: trunk supports, head supports, lateral hip supports, lumbar supports, other positioning belts

Step 9: Bring client into standing