

Power assist

# SmartDrive MX2+



permobil



Universal padded handle

Multiple ergonomic handles to quickly attach to almost any manual wheelchair

13.2 lb — Lightest rear power assist in the world

OmniWheel's patented design engages and stays engaged in any position, maintains optimal traction, and climbs 10 degree inclines.






## On-demand power assist

Pushing a manual wheelchair can put stress on shoulders and lead to serious issues in the future. SmartDrive offers a range of controls and activity-tracking technology to deliver a seamless power assist experience that may reduce repetitive stress.

Paired with the PushTracker E2 and SwitchControl, SmartDrive offers added flexibility, convenience and control options to power your push, your way.

A photograph of a paved path with scattered leaves, leading towards a dense line of green trees. In the bottom right corner, the tips of a person's feet wearing sandals are visible.

“ There were so many times I said no to something because I didn’t feel like I could. With the SmartDrive I can go all around the farmers market, go shopping, meet a friend for coffee, walk the dog and still do a ton of other things when I get home.

- Julie M., *SmartDrive Expert*


## Choose to do more

SmartDrive empowers your active lifestyle with a double tap or single touch of a button. It conserves the energy that you would normally have to use just to get around, and helps you to save it for activities that enhance your day.



Airplane travel ready



A top-down view of a person with short brown hair, wearing a light blue button-down shirt, sitting in a black wheelchair. The person's right arm is extended towards a device mounted on the wheelchair, which has a glowing blue light. The wheelchair is on a light-colored, textured floor. The background is a plain, light-colored wall.

Touch screen technology and  
bright color display makes  
the PushTracker E2 more  
accessible than ever.


8x more reliable  
Bluetooth® connection




Set up in minutes, no phone or wifi required

## PushTracker E2


Developed with connectivity in mind, the PushTracker E2 has been rigorously tested and proven to be consistently reliable and extremely durable. The PushTracker E2 makes using your SmartDrive easier than ever. With gesture-activated controls and activity-tracking technology you can start your SmartDrive with a simple tap of your hand—no pushing or reaching necessary.




Turn power assist on and off



Easily change speed settings



Review system usage



## PushTracker App

SmartDrive provides daily mobility insights for the manual wheelchair user. The PushTracker app produces and delivers activity statistics designed to promote optimal mobility health. Track distance, pushes and coast time in order to increase coast time and decrease pushes. The app can also be used to easily update and control SmartDrive settings.







SmartDrive

## SwitchControl

SwitchControl is a simple way to control your SmartDrive. Engage your SmartDrive with a push of a button for a momentary burst or press and hold to set your cruise speed for longer distances. This empowers a more diverse group of people to experience SmartDrive than ever before.






**PushTracker E2**  
Easy to read, touch  
screen settings control

**SwitchControl**  
Wired buttons with additional control  
options for your SmartDrive

**SmartDrive MX2+**  
Powers you up hills and ramps, over  
thick carpet, and long distances.

# Power your push, Your way

SmartDrive provides manual wheelchair users more control over their experience than ever before. Discover how our variety of hardware and software solutions can be combined to better support your individualized mobility goals.

	Stamina	Performance	Expert
Products available			
Best for	Hemi-propeller, pediatric, geriatric, or caregivers	Pediatric/active users (smaller, simpler wearable)	Adult/active users (advanced wearable features)
On-demand power assist	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Energy/shoulder preservation	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Easy on/off buttons	<input checked="" type="radio"/>	<input checked="" type="radio"/>	<input checked="" type="radio"/>
Hand gesture activation		<input checked="" type="radio"/>	<input checked="" type="radio"/>
Advanced Bluetooth connection			<input checked="" type="radio"/>
Easy to read, touch screen settings control			<input checked="" type="radio"/>
Performance & maintenance notifications			<input checked="" type="radio"/>





## SmartDrive Specifications

<b>Dimensions</b>	Length: 15.3" (389 mm) Width: 5.7" (141 mm) Height: 9.5" (242 mm)
<b>Product weight</b>	13.2 lb (6 kg) Drive Unit
<b>User weight</b>	30 to 331 lb (14 to 150 kg)
<b>Driving range*</b>	Up to 12.3 mi (19.8 km)
<b>Motor power</b>	250 W Brushless DC
<b>Maximum speed</b>	5.5 mph (8.8 km/h)
<b>Operating temp</b>	-13°F to 122°F (-25°C to 50°C)
<b>Chair type compatibility</b>	Rigid, Folding, Tilt-In Space, One-Arm Drive and others Chair widths: 8" (no anti-tips) -22"
<b>Drive wheel diameter</b>	22", 24", 25", 26" (501, 540, 559, 590 mm)
<b>Warranty</b>	2-Year Limited Warranty
<b>Bluetooth LE</b>	FCC ID QOQBT113 IC ID 5123A-BGTBLE113
<b>HCPCS</b>	E0986

\*Performance measurements will vary based on user and chair characteristics, driving and battery conditions.

The SmartDrive has been tested and conforms to all applicable requirements of ANSI/RESNA Standards for Wheelchairs - Volume 1 and Volume 2 and EN 12184.

NOTE: A wheelchair's specifications provided by its manufacturer may be slightly affected by the addition of the SmartDrive.

Please refer to your owner's manual for more detailed operation and care instructions.

## PushTracker E2 Specifications

<b>Dimensions</b>	1.85" x 2.06" x 0.5" (46.9 x 52.2 x 12.9 mm)
<b>Operating system</b>	Wear OS by Google™
<b>Nominal operating voltage</b>	3.9 V
<b>Electric charge</b>	415 mA-h
<b>Wireless operating band</b>	Bluetooth / BLE 2.4 GHz (2.402 to 2.480 GHz) Wi-Fi 802.11 b/g/n 2.4 GHz (2.412 to 2.472 GHz) GPS / GLONASS / Beidou 1575 MHz (1.559 to 1.610 GHz)
<b>Bluetooth LE</b>	FCC ID 2AP42-WG12026 IC ID 24006-WG12026
<b>Water resistance</b>	5 ATM (Swim-ready, up to 50 m)
<b>Band width</b>	0.87" (22 mm)
<b>Band material</b>	Silicone (interchangeable), Latex-Free and completely Biocompatible
<b>Case material</b>	Polycarbonate
<b>Display</b>	AMOLED (400 x 400 pixel)

\*All figures shown are for illustration purposes only and actual products may vary. Wear OS by Google, iPhone, Android, Bluetooth and Google Play are trademarks or registered trademarks of their respective owners.

Since U.S. Medicare coding and coverage criteria is subject to change, the provider should always confirm the HCPCS code and coverage criteria as part of the client assessment process. ©2019 Permobil. Permobil has a policy of continual product improvement and reserves the right to amend this document.