



PDG Mobility™

Comparison matrix

permobil

PDG
MOBILITY



	Fuze™ T50	Fuze T50n (NEW)	Fuze T50 Junior (peds)	Stellar GL	Stellar GLT	Stellar Impact	Stellar HD	Stellar Leap	Bentley™	Bentley LT	Bentley LT-R (NEW)	Eclipse
Dynamic tilt range	Up to 50°	Up to 50°	Up to 50°	-5° to 45°	Up to 45°	Up to 45°	Up to 45° (25° @ 600 lbs)	-30° to 20°	Up to 20°	Up to 20°	Up to 20°	None
Seat width	13" - 24"	16" to 24"	12" - 18"	14" - 32"	14" - 22"	14" - 22"	22" - 32"	14" - 22"	14" - 32"	14" - 22"	14" - 22"	20" - 48"
Seat depth	15" - 22"	15" - 22"	13" - 18"	16" - 22"	16" - 22"	17" - 22"	18" - 22"	16" - 22"	16" - 22"	16" - 22"	16" - 22"	16" - 24"
Seat to floor height	13" - 20"	16" - 20"	13" - 20"	13" - 20"	13" - 20"	14" - 21"	16" - 18"	13" - 18"	13" - 20"	13" - 20"	12" - 19"	14" - 20"
Overall width	+7.5" - 8.5"	+4.0"	+7.5" - 8.5"	+9.5 - 10.5"	+9.5 - 10.5"	+9.5 - 10.5"	+11"	+9.5 - 10.5"	+9.5" - 10.5"	+9.5" - 10.5"	+10.5"	+7.5"
Transport weight	44 lbs	44 lbs	42 lbs	51 lbs	34 lbs	64 lbs	62 lbs	58 lbs	51 lbs	33 lbs	34 lbs	44 lbs
User weight capacity	400 lbs	400 lbs	250 lbs	450 lbs	250 lbs	250 lbs	600 lbs	250 lbs	450 lbs	250 lbs	250 lbs	1000 lbs
Frame material	Aluminum	Aluminum	Aluminum	Steel	Aluminum	Steel	Steel	Steel	Steel	Aluminum	Aluminum	Steel
Power tilt option	Yes, up to 350 lbs	Yes, up to 350 lbs	No	Yes, up to 250 lbs	Yes	No	No	Yes, anterior/posterior	Yes, up to 250 lbs	Yes	Yes	No
Take-Apart	Yes	Yes	Yes	No	No	No	No	No	No	No	No	No
Hand or foot propulsion	Hand	No	Hand	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot	Hand or Foot
Self tilt	No	No	No	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Transit tiedowns	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes	Yes, up to 300 lbs	Yes	Yes	No	Yes	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes, up to 300 lbs	Yes, up to 600 lbs

Take-Apart


The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV. Available on Fuze™ T50, Fuze T50n and Fuze JR.



Seating and Positioning

Permobil provides a wide range of seating and positioning products to address the unique needs of individuals using PDG wheelchairs.

Comfort and ROHO® offer a variety of cushions, back supports and accessories that are designed for pressure redistribution, stability, and comfort that can be tailored to meet the individual needs of users.





Our innovative chair designs deliver an enhanced quality of life for both users and caregivers. Long product life, adjustability, and durability create value for all stakeholders, and our expertise in positioning wheelchairs ensures an excellent fit every time, no matter how basic or complex.