



RE: SmartDrive Firmware v3.2

04/23/2024

On April 30th, the first of several planned SmartDrive firmware/software updates of 2024 will be released. This update will be available via the PushTracker Mobile App or the SmartDrive MX2+ Wearable Apps. SmartDrive Firmware v3.2 will expand the zero position for the SpeedControl Dial, resulting in a larger window for a successful stop of the unit. As a result of this update, users may need to rotate the dial forward several additional clicks when initiating movement.

If you would like to learn how to successfully complete a firmware update, please click [here](#). Stay tuned for more product improvements coming later this year.