


TiLite manual wheelchairs | Creative Mobility Division

Aesthetic Options



How to order

Catalog

1. Identify CMDs you wish to order and fill out. Please confirm compatibility with the chair you're ordering.
2. Fill out the catalog page corresponding to that CMD with the required information for us to complete your request.
 - EZ-Ti ID, PO#, or Quote# if available.
 - Measurements, selections, and special instructions if necessary.
3. Email the filled out page(s) to chairs.tilite@permobil.com and include Order Form if applicable.

Order Forms

1. Identify CMDs you wish to order and fill out. Please confirm compatibility with the chair you're ordering.
2. In the "Notes" section on the last page of the order form, include the following information:
 - CMD number(s) being requested.
 - Measurements, selections, and special instructions if necessary.
3. Submit Order Form via email (chairs.tilite@permobil.com) or fax (509) 586-2413.

EZ-Ti

1. Identify CMDs you wish to order and fill out. Please confirm compatibility with the chair you're ordering.
2. In the "CMD" section on the Summary page of EZ-Ti, include the following information:
 - CMD number(s) being requested.
 - Measurements, selections, and special instructions if necessary.
3. Submit quote request via EZ-Ti.

Please note, custom requests may impact quote turnaround time and chair build time. Pricing for these custom options may vary based on the amount of added labor and material costs.

Unique frame/component finish including color, hydrographic and one-off tattoo design CMD #36



EZ-Ti ID _____

PO# _____

Clinical rationale

This modification is required for individuals who are requiring or requesting a personalized, unique frame color for their manual wheelchair.

Available on:

TiLite TR, TiLite TRA, TiLite Aero T, TiLite TX1, TiLite ZR, TiLite ZRA, TiLite Aero Z, TiLite 2GX, TiLite Aero X

Notes:

- Transparent film requests will require a base color call out from the customer.
- Colors are general callouts, specific color matching can be done with additional charge.
- May not work with all chair configurations.
- Note that hydrographic film pattern will stretch to varying degrees and the finish is subject to having small flaws throughout.

Special Instructions



⚠ NOTE: Permobil's warranty does not cover Creative Mobility Division frames and frame components, unless Permobil determines that the frame or frame component was defectively manufactured.

Armed forces patch(es)

CMD #74

EZ-Ti ID _____

PO# _____

Clinical rationale

Feature offered to US VA service personnel, in honor of their service, in recognition of their sacrifice.

Available on:

TiLite TR, TiLite TRA, TiLite ZR, TiLite ZRA, TiLite Aero T, TiLite Aero Z, TiLite Aero X, TiLite 2GX, TiLite TWIST.

Notes:


- Armed forces patches will be centered on back upholstery.
- Branch logo patch will be underneath Combat Wounded patch if applicable.

Required information

(Insert required information for CMD-E to complete this request)

Branch: _____

Combat Wounded: ☐ Yes OR ☐ No



Special Instructions



⚠ NOTE: Permobil's warranty does not cover Creative Mobility Division frames and frame components, unless Permobil determines that the frame or frame component was defectively manufactured.

Full polished frame CMD #87



EZ-Ti ID _____

PO# _____

Clinical rationale

Modification provides a polished appearance to titanium frame, to remove risk of paint wear/tear scratches occurring during daily routines. Provides easy buffing option during care/maintenance routine.

Available on:

Titanium frames only

Notes:

- If a Rigid wheelchair is ordered with a Titanium Footrest and Titanium Backrest, these items will be polished as well.
- If a TX1 or 2GX or Rigid 2GX, the hangers will be polished as well as the main frames as applicable.

Special Instructions



⚠ NOTE: Permobil's warranty does not cover Creative Mobility Division frames and frame components, unless Permobil determines that the frame or frame component was defectively manufactured.

Raw finished titanium frame

CMD #92

EZ-Ti ID _____

PO# _____

Clinical rationale

Modification requested by consumers who are at higher risk of wear/tear frame scratches related to ADL or MRADL routines either due to frame contact in narrow environments, or during vehicle wheelchair loading/unloading routine.

Available on:

TiLite TR, TiLite TRA, TiLite ZR, TiLite ZRA

Notes:

- Only available on Titanium frames.
- Titanium frame components will be left with surface anomalies, scuffs, weld halos and other inconsistencies that occur during the manufacturing process.



Special Instructions



⚠ NOTE: Permobil's warranty does not cover Creative Mobility Division frames and frame components, unless Permobil determines that the frame or frame component was defectively manufactured.

Attendant wheel locks for Pilot CMD #93

EZ-Ti ID _____

PO# _____

Clinical rationale

Modification indicated for safety, to allow caregivers consistent access to wheel locks from the rear of the wheelchair in tight spaces, when front access is limited (transportation vehicles).

Available on:

TiLite Pilot

Notes:

- This adaptation requires an extra-long bar on a Unilock wheel lock mounted with a 1" offset wheel lock clamp.



Special Instructions



⚠ NOTE: Permobil's warranty does not cover Creative Mobility Division frames and frame components, unless Permobil determines that the frame or frame component was defectively manufactured.

TiLite manual wheelchairs

Creative Mobility Division

2701 West Court St, Pasco, WA 99301
Phone +1 509 586 6117 | 800 545 2266



permobil