

Power wheelchairs

Power standing



permobil

Stand without compromise

Standing is an effective way to promote greater independence and reduce the risk for secondary complications associated with long-term sitting in wheelchairs. Permobil power standing wheelchairs empower you to utilise standing for activities of daily living, including the ability to drive while standing, in social, work and educational settings to positively impact mental health and quality of life. With fully programmable standing sequences that allow your position to be optimized using supine, semi-reclined, or sit-to-stand techniques.



M Corpus® VS

Midwheel Standing

The M Corpus VS combines a compact base and suspension that provides stability, traction, and ride comfort for a best-in-class power standing experience with a mid-wheel base. Designed with the user in mind, the M Corpus VS allows easy adjustment of the standing sequence & positioning to your needs in real time.

Standing &
0-300 mm Active Height

Fully customisable
standing sequences

Swing-away or
single-post knee
support options

— Option available on both
F5 Corpus VS & M Corpus VS

— Option available only on
specific chair

Functional
drive speeds
in standing



F5 Corpus® VS

Frontwheel Standing

The F5 Corpus VS has a fully independent suspension on all wheels offering comfort, traction and stability, especially beneficial for active outdoor driving. Easy adjustment of the standing sequences allows your position to be optimised while the light, swing-away knee supports offer improved comfort and flexible on-chair storage for greater independence.

Adjustable chest support

Full range of power tilt, recline, & legrest elevation

Standing &
0-350 mm Active Height

Low, functional seat-to-floor heights

Articulating power footplate



Designed for independence

With a lightweight design, our swing-away knee supports can be moved and handled independently using the easy to reach release mechanism even with limited hand function. Easily fine tune the pad placement as needed with the integrated BOA® dial* for optimal fit and comfort to your individual needs.

Available in two pad sizes, choose the option that best matches you.

*BOA® Fit System delivers micro-adjustable, precision fit in a simple to use dial



Standing supports

face-to face interactions and enables improved health outcomes.

The impact that power standing wheelchairs can have include multiple health-related benefits, as well as what activities you can participate in and how independent you are.


Identified benefits for using a power standing wheelchair supported with high quality evidence include:

- Pressure redistribution
- Bladder, bowel, digestive functions
- Range of motion and muscle length
- Activities of daily living, personal/medical care, and leisure activities
- Functional reach
- Communication
- Participation in social, work, and educational settings
- Independence

Reference: Nooijen C, Ostendorf P, Bernstein J (2023) Permobil White Paper: A summary of the evidence on power standing wheelchairs



No additional hardware required



Programming simplified

Programming changes are faster and more intuitive with QuickConfig, the new wireless programming app for Permobil power wheelchairs equipped with Power Platform. Designed for service technicians, clinicians, and providers that offers a simpler way to configure parameters and an intuitive interface to set up F5 Corpus VS and M Corpus VS with the user's preferred driving and seating characteristics and standing sequence.

Scan to
learn more



Available to download from
the App store and Google Play



App Store is a trademark of Apple Inc. Google Play is a trademark of Google LLC.

permobil

permobil.com

GR3903-01EN230531