

ATTENTION!

PLEASE HANDLE WITH CARE

Lock/Unlock - Use the lever(s) located on the wheelchair base, front or side, to lock/unlock wheels. **Use caution** as this will allow chair to roll freely.

Do not force or push the chair when in lock mode.

When lifting the chair, **only carry from the base**. Wheelchair is bottom heavy.

Do not place chair on it's side. Sensitive components can be damaged. Always keep chair standing upright.

Ensure that main fuse is switched off. Ask owner of chair for fuse's location.

Bring wheelchair to **cabin door**.

For questions, always consult the owner of the chair.

CONTACT AND WHEELCHAIR INFO

Contact Name: _____

Contact phone number: _____

Contact e-mail: _____

Chair model: _____ Colour(s) of chair: _____

Chair measurements: Weight: _____ Height: _____ Width: _____

Battery type: _____

Please refer to your owner's manual for additional safety tips and transportation guidelines.